

DcR1 antibody

Cat. No. GTX70466

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Human

Package
100 µl

Applications

Application Note

Not tested. The investigator should determine the optimal dilution for a specific application such as immunohistochemistry.

Calculated MW 27 kDa. ([Note](#))**Product Note** Reacts with the N terminus of the mature human DcR1 protein.

Properties

Form Liquid**Buffer** Serum**Preservative** 0.02% Sodium azide**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Purification** Serum**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).