

FBXO44 antibody

Cat. No. GTX70479

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

Western blot - Use at 1:500 to 1:1,000 dilution. However, the investigator should determine the optimal dilution for a specific application such as immunohistochemistry.

Calculated MW	30 kDa. (Note)
Product Note	Reacts with residues 103-119 (EWKVEDLSRDQRKEFPN) of the human FBXO44 protein. High cross reactivity: 94% identical to mouse and rat FBG3.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	residues 103-119 (EWKVEDLSRDQRKEFPN) of the human FBXO44 protein.
Purification	Affinity purified IgG
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx70479 WB Image

WB analysis of recombinant human FBXO44 protein using GTx70479 FBXO44 antibody.

Loading : 30 µg



For full product information, images and publications, please visit our [website](#).