

beta 3 Adrenergic Receptor antibody

Cat. No. GTX70685

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	ICC/IF, IHC-P
Reactivity	Human, Rat

References (2)

Package

25 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
IHC-P	10 - 15 µg/ml

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic 20 amino acid peptide from N-terminal extracellular domain of human ADRB3.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).


DATA IMAGES


GTX70685 IHC-P Image


IHC-P analysis of human adipose tissue using GTX70685 beta 3 Adrenergic Receptor antibody.
Antigen retrieval : Heat-induced antigen retrieval


GTX70685 IHC-P Image

IHC-P analysis of adrenal gland tissue using GTX70685 beta 3 Adrenergic Receptor antibody.


GTX70685 IHC-P Image

IHC-P analysis of human prostate tissue using GTX70685 beta 3 Adrenergic Receptor antibody.
Antigen retrieval : Heat-induced antigen retrieval


GTX70685 IHC-P Image

IHC-P analysis of adipose tissue using GTX70685 beta 3 Adrenergic Receptor antibody.


GTX70685 IHC-P Image

IHC-P analysis of human small intestine tissue using GTX70685 beta 3 Adrenergic Receptor antibody.
Antigen retrieval : Heat-induced antigen retrieval



For full product information, images and publications, please visit our [website](https://www.genetex.com).