

Glucagon like peptide 2 receptor antibody

Cat. No. GTX70765

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Application	IHC-P
Reactivity	Human

Package

25 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	10 µg/ml
Not tested in other applications.	

Calculated MW 63 kDa. ([Note](#))

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic 16 amino acid peptide from C-terminus of human GLP2R.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX70765 IHC-P Image

IHC-P analysis of human skin, melanoma tissue using GTX70765 Glucagon like peptide 2 receptor antibody.
Antigen retrieval : Heat-induced antigen retrieval




GTX70765 IHC-P Image

IHC-P analysis of brain, amygdala tissue using GTX70765 Glucagon like peptide 2 receptor antibody.



GTX70765 IHC-P Image

IHC-P analysis of human small intestine, intraepithelial neuroendocrine cells and enterocytes tissue using GTX70765 Glucagon like peptide 2 receptor antibody.
Antigen retrieval : Heat-induced antigen retrieval



GTX70765 IHC-P Image

IHC-P analysis of brain, pituitary tissue using GTX70765 Glucagon like peptide 2 receptor antibody.



For full product information, images and publications, please visit our [website](#).