

## PTPRA antibody

**Cat. No. GTX70943**

<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	IHC-P
<b>Reactivity</b>	Human, Mouse, Rabbit, Bovine, Dog, Hamster, Pig, Monkey, Bat, Horse

Package

25 µg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	10 µg/ml

Not tested in other applications.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.1% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	1 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Synthetic 17 amino acid peptide from cytoplasmic domain of human PTPRA.
<b>Purification</b>	Purified by affinity chromatography
<b>Conjugation</b>	Unconjugated

**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


## DATA IMAGES

**GTX70943 IHC-P Image**

IHC-P analysis of human heart tissue using GTX70943 PTPRA antibody.

Antigen retrieval : Heat-induced antigen retrieval

Dilution : 10 µg/ml

**GTX70943 IHC-P Image**

IHC-P analysis of human skeletal muscle tissue using GTX70943 PTPRA antibody.

Antigen retrieval : Heat-induced antigen retrieval

Dilution : 10 µg/ml



For full product information, images and publications, please visit our [website](#).