

GPCR150 antibody

Cat. No. GTX71087

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

Package

25 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	3 - 6 µg/ml

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic 20 amino acid peptide from 1st transmembrane domain of human GPR160.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX71087 IHC-P Image

IHC-P analysis of human colon, epithelium tissue using GTX71087 GPCR150 antibody.
Antigen retrieval : Heat-induced antigen retrieval



GTX71087 IHC-P Image

IHC-P analysis of brain, cortex, neurons and glia tissue using GTX71087 GPCR150 antibody.



GTX71087 IHC-P Image

IHC-P analysis of brain, glioblastoma tissue using GTX71087 GPCR150 antibody.



For full product information, images and publications, please visit our [website](#).