

GPR105 antibody

Cat. No. GTX71177

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	ICC/IF, IHC-P
Reactivity	Human

References (1)

Package

25 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	5 - 10 µg/ml
IHC-P	4 µg/ml

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic 20 amino acid peptide from C-terminus of human P2RY14 / P2Y14.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX71177 IHC-P Image


IHC-P analysis of human ovary, carcinoma tissue using GTX71177 GPR105 antibody.
Antigen retrieval : Heat-induced antigen retrieval


GTX71177 IHC-P Image


IHC-P analysis of human adrenal, pheochromocytoma tissue using GTX71177 GPR105 antibody.
Antigen retrieval : Heat-induced antigen retrieval


GTX71177 ICC/IF Image

ICC/IF analysis of HEK293 transfected with GPR105 using GTX71177 GPR105 antibody.


GTX71177 IHC-P Image

IHC-P analysis of adrenal gland tissue, pheochromocytoma using GTX71177 GPR105 antibody.


GTX71177 IHC-P Image

IHC-P analysis of breast, adenocarcinoma tissue using GTX71177 GPR105 antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).