

GPR44 antibody

Cat. No. GTX71385

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	ICC/IF, IHC-P
Reactivity	Human

Package

25 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
IHC-P	5 - 10 µg/ml

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic 33 amino acid peptide from N-terminal extracellular domain of human GPR44 / CRTH2.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX71385 IHC-P Image

IHC-P analysis of lung, small cell carcinoma tissue using GTX71385 GPR44 antibody.


GTX71385 IHC-P Image


IHC-P analysis of human nasal mucosa, eosinophils tissue using GTX71385 GPR44 antibody.

Antigen retrieval : Heat-induced antigen retrieval



GTX71385 IHC-P Image

IHC-P analysis of human lymph node, non-hodgkins lymphoma tissue using GTX71385 GPR44 antibody.

Antigen retrieval : Heat-induced antigen retrieval


GTX71385 IHC-P Image

IHC-P analysis of human tonsil tissue using GTX71385 GPR44 antibody.


GTX71385 IHC-P Image

IHC-P analysis of brain, thalamus, neuron tissue using GTX71385 GPR44 antibody.



For full product information, images and publications, please visit our [website](#).