

GPR78 antibody

Cat. No. GTX71513

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

Package

25 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	14 µg/ml

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic 17 amino acid peptide from 3rd cytoplasmic domain of human GPR78.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX71513 IHC-P Image


IHC-P analysis of human placenta, villi tissue using GTX71513 GPR78 antibody.
Antigen retrieval : Heat-induced antigen retrieval


GTX71513 IHC-P Image


IHC-P analysis of human skin, melanoma tissue using GTX71513 GPR78 antibody.
Antigen retrieval : Heat-induced antigen retrieval


GTX71513 IHC-P Image

IHC-P analysis of brain, alzheimer's disease, neurofibrillary tangles tissue using GTX71513 GPR78 antibody.


GTX71513 IHC-P Image

IHC-P analysis of brain, alzheimer's disease, senile plaque tissue using GTX71513 GPR78 antibody.


GTX71513 IHC-P Image

IHC-P analysis of brain, alzheimer's disease, neurofibrillary tangles tissue using GTX71513 GPR78 antibody.



For full product information, images and publications, please visit our [website](#).