

Somatostatin antibody [SOM-018]

Cat. No. GTX71935

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ICC/IF, IHC-P, IHC-Fr, IHC
Reactivity	Human, Mouse, Rat

References (56)

Package

500 µl

Applications

Application Note

Recommended Starting Dilutions:

For ICC/IF: Use at a dilution of 1:25-1:100

For IHC-Fr: Use at a dilution of 1:25-1:100

For IHC-P: Use at a dilution of 1:25-1:100

Not yet tested in other applications. Optimal dilutions should be determined experimentally by the researcher.

Product Note The epitope recognised by this antibody lies within the Somatostatin-14.

Properties

Form	Liquid
Buffer	PBS, 0.2% BSA
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	0.14 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Somatostatin conjugated to a protein carrier
Purification	Protein A affinity purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx71935 IHC-Fr Image

The data was published in the journal J Neurosci in 2019. [PMID: 30926750](https://pubmed.ncbi.nlm.nih.gov/30926750/)


GTx71935 IHC-Fr Image

The data was published in the journal J Neurosci in 2019. [PMID: 30926750](https://pubmed.ncbi.nlm.nih.gov/30926750/)


GTx71935 IHC Image

The data was published in the journal Hippocampus in 2017. [PMID: 27997999](https://pubmed.ncbi.nlm.nih.gov/27997999/)


GTx71935 IHC Image

The data was published in the journal Hippocampus in 2017. [PMID: 27997999](https://pubmed.ncbi.nlm.nih.gov/27997999/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).