

PGP9.5 antibody [13C4]

Cat. No. GTX72095

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human, Mouse, Rat

Package
500 µl

Applications

Application Note

For IHC-P: Use at a dilution of 1:30. It is recommended to perform heat mediated antigen retrieval before starting with IHC protocol. Optimal dilutions/concentrations should be determined by the researcher.

Properties

Form	Liquid
Buffer	PBS, 0.2% BSA
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	0.3 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	PGP 9.5 from human brain
Purification	Protein A affinity purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).