

TSP antibody [A4.1]

Cat. No. GTX72221

Host	Mouse
Clonality	Monoclonal
Isotype	IgM
Applications	WB, ICC/IF
Reactivity	Human

Package
500 µl

Applications

Application Note

For IF: Use at an assay dependent dilution. For WB: Use at a dilution of 1:100. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW 129 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.2% BSA
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Calcium-replete, native purified human TSP from the supernatant of thrombin-activated platelets
Purification	Protein A affinity purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).