

CD155 antibody [D171]

Cat. No. GTX72348

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ICC/IF, FCM, Blocking
Reactivity	Human, Monkey

References (1)

Package

500 µl

Applications

Application Note

For FACS: Use at an assay dependent dilution. For IF: Use at an assay dependent dilution. Inhibits Replication of Poliovirus in HeLa Cells. Not suitable for western blot.

Product Note Inhibits Replication of Poliovirus in HeLa Cells

Properties

Form	Liquid
Buffer	PBS, 0.2% BSA
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The epitope is : aa 35-50 from human
Purification	Protein G affinity purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).