

CDw17 antibody [B392 (HuLy-m13)]

Cat. No. GTX72455

Host	Mouse
Clonality	Monoclonal
Isotype	IgM
Applications	ICC/IF, IHC-Fr, FCM
Reactivity	Human

Package

200 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	5-20 µg/ml
IHC-Fr	1:50 to 1:100
FCM	1-1.5 µg for 10 ⁶ cells

Not tested in other applications.

Product Note This antibody reacts with lactosylceramide (cer-Glc b1->4 Gal).

Properties

Form	Liquid
Buffer	PBS, 0.2% BSA
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	0.4 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Immunoprecipitated complexes of human detergent PBL lysate and anti-b2m proteins.
Purification	Protein A affinity purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).