

Testosterone antibody

Cat. No. GTX72779

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P, IHC-Fr, IHC
Reactivity	Species independent

References (5)

Package

500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	Assay dependent
IHC-Fr	Assay dependent
IHC	Assay dependent

Note : This antibody is prediluted for Immunohistochemistry use. Prolonged fixation in buffered formalin can destroy the epitope.

Not tested in other applications.

Product Note Indirect Testosterone: 100%, Androstenedione: 0.5% and DHEA: 0.01%

Properties

Form	Liquid
Buffer	Serum
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Immunogen	Testosterone-BSA
Purification	Unpurified From polyclonal serum
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).