

## C Reactive Protein antibody

Cat. No. GTX72817

<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	IP, ELISA, RIA
<b>Reactivity</b>	Human

Package  
500 µl

## Applications

## Application Note

For ELISA, IP, and RIA: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the researcher.

**Product Note** C-reactive protein - Does not exhibit cross-reactivity against normal human serum

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	15mM Phosphate, 100mM NaCl
<b>Preservative</b>	0.1% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Purified human C-reactive protein.
<b>Purification</b>	Purified by DEAE anion exchange chromatography.
<b>Conjugation</b>	Unconjugated

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).