

C Reactive Protein antibody

Cat. No. GTX72817

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IP, ELISA, RIA
Reactivity	Human

Package
500 µl

Applications

Application Note

For ELISA, IP, and RIA: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the researcher.

Product Note C-reactive protein - Does not exhibit cross-reactivity against normal human serum

Properties

Form	Liquid
Buffer	15mM Phosphate, 100mM NaCl
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Purified human C-reactive protein.
Purification	Purified by DEAE anion exchange chromatography.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).