

LH beta antibody

Cat. No. GTX72837

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	ICC/IF, IHC
Reactivity	Human

Package
500 µl

Applications

Application Note

IHC, ICC (frozen or formalin-fixed paraffin-embedded (FFPE) tissue sections, cell smears) 1:100-1:500 using streptavidin-Biotin system or polymer system.

Product Note

This antibody stains both free and intact LH and is useful for the classification of pituitary adenomas and tumors of the pituitary.

Properties

Form	Liquid
Buffer	10mM PO ₄ , 150mM NaCl 1% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human LH beta subunit.
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).