

MUC2 antibody

Cat. No. GTX72840

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	ICC/IF, IHC
Reactivity	Human

Package
500 µl

Applications

Application Note

IHC, ICC (frozen or formalin-fixed paraffin-embedded (FFPE) tissue sections, cell smears) 1:50-1:200 using streptavidin-Biotin system or polymer system. Prediluted antibody is used directly without dilution. For FFPE tissue sections high temperature antigen retriever (boiling tissue in 10 mM citrate buffer for 10 minutes) enhances the signal and may require further dilution of antibody.

Properties

Form	Liquid
Buffer	10mM PO ₄ , 200mM NaCl, 1% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic MUC2 (23 aa) tandem repeat consensus sequence
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).