

alpha Amylase antibody

Cat. No. GTX72901

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P, ELISA, Double immunodiffusion
Reactivity	Human

Package
500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:100
ELISA	1:4,000-1:6,000
Double immunodiffusion	Assay dependent

Not tested in other applications.

Product Note

The antiserum is specific for human α-amylase found in human saliva and human pancreatic extract. No reaction with other human saliva proteins or pancreatic extract proteins is observed.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Immunogen	human salivary α-amylase.
Purification	fractionated antiserum
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).