

Avidin antibody

Cat. No. GTX72952

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA
Reactivity	Chicken

Package
1 ml

Applications

Application Note

ELISA-1,000-5,000 . Antibody Characteristics: Immunogen: Highly purified egg white avidin. Purity: The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. Specificity: The reactivity of this antibody to avidin was demonstrated by immunodiffusion techniques on agar plates. For research use only. Not intended or approved for use as diagnostic purposes.

Product Note This antibody reacts with 68 kD avidin and does not cross-react with streptavidin.

Properties

Form	Liquid
Buffer	10mM PO ₄ , 150mM NaCl
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Avidin
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).