

# Goat Anti-Rat IgG antibody, pre-adsorbed (HRP)

**Cat. No. GTX73151**

<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	WB, IHC-P, IHC-Fr, ELISA
<b>Reactivity</b>	Rat

References ( 3 )

Package

500 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
IHC-Fr	1/50
ELISA	1/1000

Not tested in other applications.

**Product Note** Cross-reactivity with mouse immunoglobulins has been removed by solid phase adsorption.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, HRP stabilizer
<b>Preservative</b>	0.01% Thimerosal
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
<b>Concentration</b>	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Rat IgG
<b>Purification</b>	Purified by affinity chromatography
<b>Conjugation</b>	Horseradish peroxidase(HRP)

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).