

Digoxigenin antibody

Cat. No. GTX73152

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA
Reactivity	Species independent

References (2)

Package

1 mg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	1:100,000-1:1,000,000

Note : as capture antibody

Not tested in other applications.

Product Note

On ELISA this antibody reacts equally with digoxigenin-BSA and digoxin-BSA; This antibody does not react with carrier protein used for conjugation.

Properties

Form	Liquid
Buffer	10mM Phosphate, 150mM NaCl
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Digoxigenin-conjugated with a carrier protein.
Purification	Purified by antigen-affinity chromatography and ion exchange chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).