

Bovine Serum (Calf) - Sterile

Cat. No. GTX73232

Application	WB, ICC/IF, ELISA, Blocking, IHC
Species	Bovine

Package
500 ml

PRODUCT

Summary	<p>Bovine Calf Serum</p> <p>U.S. sourced, ideal for those who seek an uncompromising quality</p> <p>Low hemoglobin</p> <p>Endotoxin: < 1.0 ng/ml</p> <p>Collected from animals of US origin</p>
----------------	--

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
ELISA	Assay dependent
Blocking	Assay dependent
IHC	Assay dependent

Not tested in other applications.

PROPERTIES

Form	Liquid
-------------	--------

Storage	Keep as concentrated solution. Store at -20°C or below. Avoid multiple freeze-thaw cycles. Use standard sterile procedures when handling this product to maintain sterility.
----------------	--

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
-------------	---

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).