

Bovine Plasma (Calf) - Sterile

Cat. No. GTX73238

Application	WB, ICC/IF, ELISA, Blocking, IHC
Species	Bovine

Package
500 ml

PRODUCT

Summary	Low hemoglobin. Collected from animals of US origin U.S. sourced, ideal for those who seek an uncompromising quality.
----------------	---

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

WB	Assay dependent
----	-----------------

ICC/IF	Assay dependent
--------	-----------------

ELISA	Assay dependent
-------	-----------------

Blocking	Assay dependent
----------	-----------------

IHC	Assay dependent
-----	-----------------

Not tested in other applications.

PROPERTIES

Form	Liquid
-------------	--------

Buffer	For inquiries regarding the availability of plasmas with particular anticoagulants, please contact support@genetex.com .
---------------	---

Preservative	No preservative
---------------------	-----------------

Storage	Keep as concentrated solution. Store at -20°C or below. Avoid multiple freeze-thaw cycles. Use standard sterile procedures when handling this product to maintain sterility.
----------------	--

Endotoxin	< 1.0 ng/ml
------------------	-------------

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
-------------	---

	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
--	---



For full product information, images and publications, please visit our [website](https://www.genetex.com).