

# Fetal Bovine Serum (Charcoal Stripped) - Sterile

**Cat. No. GTX73240**

**Applications** WB, ICC/IF, ELISA, Blocking, IHC

**Species** Bovine

**References** ( 2 )

**Package**

100 ml

## PRODUCT

### Summary

FBS-charcoal stripped Triple 0.1 micron sterile filtered, 9CFR TESTED, This sera has reduced levels of hormones and is utilized in receptor or estrogen studies. This treatment minimized variability in sera lot to lot and normally improves immunoassay systems. Performance tested Hemoglobin level: <20 mg/ dl Adventitious virus tested Ideal for those who seek an uncompromising quality FBS at a reasonable price.

### Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
ELISA	Assay dependent
Blocking	Assay dependent
IHC	Assay dependent

Not tested in other applications.

## Properties

**Form** Liquid

**Storage** Keep as concentrated solution. Store at -20°C or below. Avoid multiple freeze-thaw cycles. Use standard sterile procedures when handling this product to maintain sterility.

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).