

Eosin Y Counterstain (Ready-to-use)

Cat. No. GTX73303

Application IHC

Package
250 ml

PRODUCT

Summary Eosin Y Counterstain - 250 ml

APPLICATION

Application Note

Please see protocol for details.

PROPERTIES

Form Liquid

Storage Store at room temperature.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).