

# Light Green counterstain (Ready-to-use)

**Cat. No. GTX73306**

**Applications**

IHC

References ( 1 )

Package

50 ml

## PRODUCT

### Summary

Light Green Counterstain

### Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

### Suggested dilution

### Recommended dilution

IHC

Assay dependent

Not tested in other applications.

## Properties

### Storage

Store at room temperature.

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).