

Light Green counterstain (Ready-to-use)

Cat. No. GTX73306

Applications

IHC

References (1)

Package

50 ml

PRODUCT

Summary

Light Green Counterstain

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution

Recommended dilution

IHC

Assay dependent

Not tested in other applications.

Properties

Storage

Store at room temperature.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).