

Tissue conditioner

Cat. No. GTX73312

Application IHC

Package
15 ml

PRODUCT

Summary Tissue Immunostain Conditioner

APPLICATION

Application Note

This product is designed to unmask immunoreactive sites altered by fixation and/or the embedding process. Please see application protocol.

PROPERTIES

Storage Keep as concentrated solution, aliquot and store at 4°C. DO NOT FREEZE.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).