

Brij 35 (Polyoxyethylene Laural Ether)

Cat. No. GTX73318

Application WB, IHC-P

Package
50 ml

PRODUCT

Summary Brig 35 (Polyoxyethylene Laural Ether) - 50ml

APPLICATION

Application Note

For diluting antibodies use 1ml Brig 35 solution added to 400 ml of appropriate buffer. For washing buffer use 5 ml Brig 35 solution added to 95 ml of appropriate buffer. Optimal dilutions/concentrations should be determined by the individual lab.

PROPERTIES

Form Liquid

Buffer This detergent is supplied as 15% in water

Preservative 0.15 % Proclin 300

Storage Store at room temperature.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).