

10X ImmunoStain wash buffer

Cat No. GTX73344

Application IHC

Reference (1)
 Package
 1 liter

PRODUCT

Summary 10X ImmunoStain wash buffer (pH 7.5)

APPLICATION

Application Note

No pH adjustment or additives are required when using this formulation. Simply dilute the concentrated 10-fold buffer with distilled or deionized water. A total of 10 liters of working rinsing solution can be prepared with 1 liter of the concentrated reagent. This buffer is also the reagent of choice for rinsing tissue sections that are subjected to immunocytochemical staining using conventional horizontal processing techniques.

PROPERTIES

Storage Store at room temperature.

Note For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).