

C3 / C3c antibody (FITC)

Cat. No. GTX73472

Host	Rabbit	
Clonality	Polyclonal	
Isotype	IgG	
Applications	ICC/IF, IHC-P	
Reactivity	Human, Mouse, Rat, Goat, Sheep, Bovine, Cat, Dog, Guinea pig, Pig, Mink	

References (1)
Package
1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution	
ICC/IF	Assay dependent	
IHC-P	1:75-1:100	

Note: Antigen retrieval: Proteinase K 30 min at 37°C

Not tested in other applications.

Product NoteThis antibody reacts with human C3c complement and with the C3c part of C3 and C3b.

Properties	
Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	C3c Complement isolated from complement-activated human serum.
Purification	Purified IgG
Conjugation	Fluorescein isothiocyanate (FITC) Wavelength
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 07 Page 1 of 1