

## Human Papillomavirus type 16 L1 antibody [CamVir 1]

Cat. No. GTX73652

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ICC/IF, IHC-P, IP
Reactivity	Human papillomavirus

References ( 2 )

Package

1 ml

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-P	1:25-1:50
IP	Assay dependent

**Note : Recommendations for antigen retrieval : EDTA buffer pH 8.0**

Not tested in other applications.

## Properties

Form	Liquid
Buffer	Ascites
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Immunogen	BALB/C mice were immunized with recombinant HPV-16 L1 protein.
Purification	Unpurified
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX73652 IHC-P Image

IHC-P analysis of human HPV type 16 infected cervix tissue using GTX73652 Human Papillomavirus type 16 L1 antibody [CamVir 1].



For full product information, images and publications, please visit our [website](#).