

## CD79a antibody [HM57]

Cat. No. GTX74022

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Application	IHC-P, IHC-Fr, FACS
Reactivity	Human, Mouse, Rat, Rabbit, Bovine, Chicken, Guinea pig, Pig, Horse, Opossum

Reference ( 1 )  
Package  
100 µg

## APPLICATION

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FACS	Assay dependent

Not tested in other applications.

Calculated MW 25 kDa. ( [Note](#) )

## PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide corresponding to amino acids 202-216 of human CD79a
Purification	Protein A purified
Conjugation	Unconjugated

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX74022 FACS Image

FACS (intracellular staining) analysis of human peripheral blood using GTX74022 CD79a antibody [HM57].



For full product information, images and publications, please visit our [website](#).

Date 2024 / 05 / 20 Page 2 of 2