

CXCL1 / GRO alpha antibody

Cat. No. GTX74108

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP
Reactivity	Mouse, Rat

References (1)

Package

50 µg

Applications

Application Note

For WB: Tested using transfected cell lines only. Use at a dilution of 1:2,000-1:5,000. For IP: Use at a dilution of 1:300-1:800. Not suitable for use in IHC. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW 10 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment - amino acids 1-72 of rat KC
Purification	Protein A affinity purified From polyclonal serum
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX74108 WB Image

Rat KC detected by WB in the supernatant of KC-transfected cells.



For full product information, images and publications, please visit our [website](#).