

Transferrin antibody

Cat. No. GTX74139

Host	Sheep
Clonality	Polyclonal
Isotype	IgG
Applications	Immunodiffusion
Reactivity	Human, Mouse, Rat, Rabbit, Cat, Guinea pig

Package 1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Immunodiffusion	Neat

Not tested in other applications.

Properties	
Form	Liquid
Buffer	Glycine buffered saline, 0.1% EACA, 0.01% Benzamidine, 1mM EDTA
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Human holo-transferrin.
Purification	Purified by ion exchange chromatography
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 06 Page 1 of 1