

Albumin antibody

Cat. No. GTX74141

Host	Sheep
Clonality	Polyclonal
Isotype	IgG
Application	Immunodiffusion
Reactivity	Human, Mouse, Rat, Rabbit, Bovine, Cat, Dog, Chicken, Guinea pig, Pig, Horse

Package
1 ml

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Immunodiffusion	Assay dependent
Not tested in other applications.	

Calculated MW 69 kDa. ([Note](#))

PROPERTIES

Form	Liquid
Buffer	Glycine buffered saline, 0.1%EACA0.01% Benzamidine, 1mM EDTA
Preservative	0.099% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Human serum albumin purified from Cohn fraction V.
Purification	Purified by ion exchange chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).