

ICAD antibody

Cat. No. GTX74223

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, ELISA
Reactivity	Human

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1µg/ml
ICC/IF	5µg/ml
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	37 kDa. (Note)
Product Note	This antibody recognizes non-cleaved (45kDa) and cleaved DFF45.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide corresponding to aa 313-331 at the C-terminus of human DFF45 (accession no. NP_004392).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX74223 WB Image

Western blot analysis of whole cell lysates from HeLa (H), Jurkat (J), A431 (A) and K562 (K) cell lines probed with Rabbit anti Human DFF (GTX74223)




GTX74223 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX74223 ICAD antibody.

Green : Primary antibody

Blue : DAPI

Dilution : 20 µg/ml



GTX74223 ICC/IF Image

Immunocytochemical staining of HeLa cells with Rabbit anti Human DFF (GTX74223)



For full product information, images and publications, please visit our [website](https://www.genetex.com).