

ICAD antibody

Cat. No. GTX74224

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP
Reactivity	Human

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	2µg/ml
IP	Assay dependent

Not tested in other applications.

Calculated MW 37 kDa. ([Note](#))

Product Note This antibody recognizes non-cleaved (45kDa) and cleaved DFF.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide corresponding to aa 2-21 at the N-terminus of human DFF45 (accession no. NP_004392).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX74224 WB Image

WB analysis of HeLa (H), K562 (K), Jurkat (J), and Raji (R) cell lysates using GTX74224 ICAD antibody.

Dilution : 2µg/ml



For full product information, images and publications, please visit our [website](#).