

Caspase 3 antibody

Cat. No. GTX74264

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

References (1) Package 50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

*Optimal dilutions/concentrations should be determined by the researcher.		
Suggested dilution	Recommended dilution	
WB	1/1,000	
IHC-P	1/100-1/500	
Note: This product requires antigen retrieval using heat treatment prior to staining of paraffin section.		
Not tested in other applications.		
Calculated MW	32 kDa. (<u>Note</u>)	
Product Note	This antibody recognizes both intact Caspase 3 and the 17kD fragment, but not the 11kD fragment.	
Properties		
Form	Liquid	
Buffer	Serum	
Preservative	0.09% Sodium azide	
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.	
Immunogen	Recombinant human caspase-3.	
Purification	Serum	
Conjugation	Unconjugated	
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.	
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to	



For full product information, images and publications, please visit our <u>website</u>.

determine the structure or sequence of the product.

Date 2025 / 11 / 07 Page 1 of 1