

CD40 antibody [3/23]

Cat. No. GTX74689

Host	Rat
Clonality	Monoclonal
Isotype	lgG2a
Applications	IHC-Fr, FCM, ELISA
Reactivity	Mouse

Package 100 μg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution	
IHC-Fr	Assay dependent	
FCM	1/10-1/50	
ELISA	Assay dependent	
Note: the 10st of the competed westing distinct to lebel 106 cells in 100st		

Note : Use $10\mu l$ of the suggested working dilution to label 10^6 cells in $100\mu l$.

Not tested in other applications.


Properties	
Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Extracellular Domain of Mouse CD40 and the Fc portion of Human IgG1.
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 04 Page 1 of 2

DATA IMAGES



GTX74689 FCM Image

FACS analysis of mouse splenocytes using GTX74689 CD40 antibody [3/23].



For full product information, images and publications, please visit our website.

Date 2025 / 11 / 04 Page 2 of 2