

Mouse IgG2a isotype control (FITC)

Cat. No. GTX74757

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	FCM

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	Assay dependent

Not tested in other applications.

Product Note

This antibody is suitable for use as a negative control to assess non-specific binding of mouse IgG2a antibodies to rat cells. And It is not suitable for use as a negative control on human.

Properties

Form Liquid

Buffer PBS, 1% BSA

Preservative 0.09% Sodium azide

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light.

Concentration 0.1 mg/ml (Please refer to the vial label for the specific concentration.)

Purification Protein A purified

Conjugation Fluorescein isothiocyanate (FITC) [Wavelength](#)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).