

## CD62L antibody [MEL-14]

Cat. No. GTX74799

Host	Rat
Clonality	Monoclonal
Isotype	IgG2a
Applications	FCM, IP, IHC
Reactivity	Mouse

Package  
250 µg

## Applications

## Application Note

Optimal dilutions/concentrations should be determined by the researcher.

Product Note	Mouse CD62L
--------------	-------------

## Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Mouse B cell lymphoma, 38C-14
Purification	Protein G purified
Conjugation	Unconjugated

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTx74799 FCM Image

Staining of mouse spleen with Rat anti Mouse CD62L: RPE (GTx74800)



For full product information, images and publications, please visit our [website](#).