

## CD161 antibody [10-78] (PE)

Cat. No. GTX75001

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	FCM
Reactivity	Rat

**Package**  
100 test

## Applications

## Application Note

For FACS: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the researcher.

## Properties

Form	Liquid
Buffer	PBS, BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified splenic NK cells from the LEW rat strain.
Purification	Purified IgG
Conjugation	Phycoerythrin (PE) <a href="#">Wavelength</a>


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX75001 FCM Image**

Staining of rat peripheral blood lymphocytes with Mouse anti RatT CD161(GTX74998):Alexa Fluor® 647



For full product information, images and publications, please visit our [website](#).

Date 2026 / 02 / 02 Page 2 of 2