

# CD19 antibody [6D5] (PE-Cy5)

**Cat. No. GTX75014**

<b>Host</b>	Rat
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2a
<b>Applications</b>	FCM
<b>Reactivity</b>	Mouse

**Package**  
50 test

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	Neat-1/5

**Note :** Use 10µl of the suggested working dilution to label 10<sup>6</sup> cells or 100µl whole blood.

Not tested in other applications.

**Product Note** Clone 6D5 has been shown to recognize the same epitope as clone ID3 in cross-competition assays.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, Sucrose
<b>Preservative</b>	0.1% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Human K562 cell line transfected with murine CD19.
<b>Purification</b>	Purified IgG
<b>Conjugation</b>	Phycoerythrin-Cyanine5 (PE-Cy5) <a href="#">Wavelength</a>

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).