

## hCG alpha antibody [INN-hCG-5]

Cat. No. GTX75049

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P, IHC-Fr, RIA
Reactivity	Human

References ( 1 )

Package

100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	Assay dependent
IHC-Fr	Assay dependent
RIA	Assay dependent

Not tested in other applications.

## Product Note

Clone INN-h-CG-5 recognizes an epitope on the alpha subunits of all human glycoprotein hormones. Indirect Immunofluorescence analysis on paraffin sections of human pituitaries reveals clear staining of glycoprotein hormone producing cells. Reacts with human LH, FSH, TSH, and hCG.

## Properties

Form	Liquid
Buffer	Ascites
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).