

## Poly (Adp-Ribose) antibody [10H]

**Cat. No. GTX75054**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG3
<b>Applications</b>	WB, IHC-P, IHC-Fr, IP
<b>Reactivity</b>	Species independent

References ( 8 )

Package

50 µl

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500
IHC-P	Assay dependent
IHC-Fr	Assay dependent
IP	Assay dependent

**Note : This product requires protein digestion pre-treatment of paraffin sections e.g. trypsin or pronase.**

Not tested in other applications.

**Product Note** The antibody does not cross-react with RNA, DNA or monomers of ADP-ribose.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, 50% Glycerol
<b>Preservative</b>	No preservatives
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Purified poly (ADP-ribose) polymer, 10-50 unit chain length
<b>Purification</b>	Purified IgG
<b>Conjugation</b>	Unconjugated



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note**

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).