

Poly (Adp-Ribose) antibody [10H]

Cat. No. GTX75054

Host	Mouse
Clonality	Monoclonal
Isotype	IgG3
Applications	WB, IHC-P, IHC-Fr, IP
Reactivity	Species independent

References (8)

Package

50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500
IHC-P	Assay dependent
IHC-Fr	Assay dependent
IP	Assay dependent

Note : This product requires protein digestion pre-treatment of paraffin sections e.g. trypsin or pronase.

Not tested in other applications.

Product Note The antibody does not cross-react with RNA, DNA or monomers of ADP-ribose.

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified poly (ADP-ribose) polymer, 10-50 unit chain length
Purification	Purified IgG
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).