

# CD71 antibody [OX-26]

**Cat. No. GTX75149**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2a
<b>Applications</b>	WB, FCM
<b>Reactivity</b>	Mouse, Rat

**Package**  
500 µg

## Applications

### Application Note

Flow cytometry, WB; Flow Cytometry: Use 10ul of the suggested working dilution to label 10<sup>6</sup> cells in 100ul.

**Product Note**      rat CD71

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	No preservatives
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	1 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	PHA activated rat lymphocytes
<b>Purification</b>	Purified IgG
<b>Conjugation</b>	Unconjugated


### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX75149 FCM Image**

Staining of mouse spleen cells with Mouse anti Rat CD71(GTX75149):Low Endotoxin


**GTX75149 FCM Image**

Staining of mouse spleen with Mouse anti Rat CD71(GTX75149 ):Biotin


**GTX75149 FCM Image**

Staining of rat spleen lymphocytes with Mouse anti Rat CD71 (GTX75149)



For full product information, images and publications, please visit our [website](https://www.genetex.com).